WS 316

Faeko

Sampler

Housing

Complete housing of stainless steel **including** roof, robust and durable

Simple upgrade to Monitoring Station

XY-distributor

Over 10 years of successful operation

Direct dosing system, **no cleaning** of distribution plate, no cross-contamination

Standard up to **36 x 1 ltr.** bottles (optionally 48 x 1 ltr. bottles)


Numerous bottle combination are realisable (7 bottle combinations are pre-installed)

Different bottle sizes composite and reference samples can be bottles at the same time

Client specific bottles can be used without difficulty

Discharge and **rinsing** allow the cleaning of the complete sampling line

Dosing



Recalibration of sample volume after cleaning is **not necessary**

Dosing vessel situated in **protected and thermostatted** sample chamber,
thus **additional** panel door and
protective heater **not necessary**

Hose inlet from the **left** or **right** optionally from the **bottom**

Closing pressure of pinch-valve is **adjustable**

Reliable sampling and accurate sample volumes with vacuum pressure system

Parallel programs and program linkage are standard





Water Sam® A versatile sampler according to ISO 5667. Including a built-in cleaning system for washing out dosing unit and hose. This makes the WS 316 Faeko well suited to difficult situations where the sample medium might contain particles that would clog a normal sampler.

MicroSam The control meets all your sampling requirements such as time, volume and event proportional sampling. A manual sample is possible at any time without interrupting the program.

Up to 9 programs can individually be set by the customer. They can be interlinked or run simultaneously.

The integrated control unit protects your equipment and logs important data such as pump run time, run time of the cooling unit, failure messages etc.

A GSM modem for remote control can be connected. An SMS error message can then also be sent to a defined mobile phone.

Examples:







1 x 25.0 litres

4x 15.0 litres

15 x 2.9 litres





24x 1.0 litre

35x 1.0 litres

With the special air circulation system inside the storage compartment the samples are kept at optimum temperature conditions of +4°C (at ambient temperatures from -25°C up to +42°C).

Electrical and mechanical components are situated in the upper dry part of the sampler.

The housing can be disassembled and recycled.

Service

Water Sam® samplers are very easy to maintain. Our customer service is always ready to help you fast and reliably.

Manufactured by:



WaterSam GmbH & Co. KG Felix-Wankel-Straße 70 72108 Rottenburg am Neckar Deutschland

Tel: + 49 (0)7472 9882-0 Fax: + 49 (0)7472 9882-22 E-Mail: info@watersam.de www.watersam.com

Technical Data - WS 316

According to ISO 5667

According to ISO 5667	
Sampling Method	Vacuum dosing system (Peristaltic pump upon request)
Material with fluid contact	Borosilicate glass, PE, PVC, stainless steel AISI 304, AISI 316 Ti, silicone
Rinsing	With clean— or process water, incl. pressure regulator and valves
Dosing volume	20 – 350 ml (multiple dosing is possible)
Suction hose	di 12 mm, PVC
Suction height	7.5 m, pump capacity 14.5 l/min without back pressure
Pump	230 V AC; -0,8 to 1 bar; no collector
Distribution	X-Y distributor, directly into the bottles
Distribution hose	di 15 mm, silicone
Housing	Stainless steel (AISI 304) Option: AISI 316 Ti or powder-coated (RAL colors) 40mm insulation, cold bridge free
Cooling unit	Compressor cooling 230 V AC; 160 W; refrigerant R 134a CFC-free; automatic defrost
Heating unit	Electric heating in stainless steel case 230 V; 350 W with excess temperature protection 70°C
Thermostat	PT100
Control	Microprocessor controlled, 4 button tactile back lighted keypad with 4 x 20 col. LCD display, Real time clock, 5 year RAM battery back-up, overload protected outputs
Software	Menu driven, 9 possible programs (simultaneously)
	Memories for errors, events and operation status
	Preset programs, easy to modify
	Preset and free distributor settings
	Combination and interlinking programs as standard features
Inputs	0/4 00 4
Flow signals : Event signals :	0/4-20mA or. potential free contact Potential free contact
Outputs standard basic model	4 outputs all in all program is running; error message 2 outputs free
Interface	RS 232 or optional RS 485
Electrical Data	Mains supply 230 V / 50 Hz Power input approx. 750 VA max. Main fuse 16 A
Dimensions	Height 1.290 mm Width 655 mm (including fastening bar =720mm) Depth 770 mm weight 90-125 kg (depending on equipment)

Supplied by: